



Plug n' Play Diagnostic Hysteroscopy Set



Adler Instruments announces a new era in hysteroscopy with the unveiling of our new Hysteroscopy Set. Utilizing our 2.7mm hysteroscope, our simplified endoscopic camera and LED light source, visualization of the entire uterine cavity can be seen with brilliant clarity. The integrated continuous-flow single sheath system simplifies usage even more. The laterally angled in-and-outflow valves guarantee optimal irrigation as well as ergonomic comfort.

We also offer: Refurbished office endoscopy systems from : Karl Storz, Olympus, Welch Allen, Richard Wolf.

- **Smaller optics and sheath design allow for easier insertion, less dilatation, and greater patient comfort.**



- **Outstandingly suitable for diagnostic and operative procedures.**



700-063



Set Includes:

- Continuous Flow Sheath
- 2.7mm Hysteroscope
- 1 Chip CCD Camera
- Coupler
- LED Light Source
- Fiberoptic Light Cable
- Choice of Instrument:
 - Alligator Forcep
 - Scissors
 - Oval Punch



250-807



214-505-3896 • 888-867-7118 Fax

adlerinstruments.com