

Small Diameter Endoscope Specialist




Certified Factory Trained Technicians,
Comprehensive Repair, Restoration and Sales
Repair Exchange Program
Loaner Scopes Available

OLYMPUS® PENTAX **STORZ**
KARL STORZ – ENDOSKOP

Gold Fiducial Marker Needles

Square Hub Marker Kit



- Electropolished for sharpness
- Siliconized for easy insertion
- Mirror polished I.D. for a smooth, easy seed deployment
- Echo-enhanced tipped for ultrasound visualization
- Centimeter markings for depth placement
- Rounded stylet for smooth gold marker delivery
- Foil marker to indicate bevel tip orientation

BARD™ Style

Magnum Style & Biopty/Promag Prostate Biopsy Needles



Pro Mag/Biopty Style
UMRPM-1820 18g x 20cm

- Centermeter markings for depth placement
- Electropolished for sharpness and smooth entry
- 19mm or 17mm cutting notch for tissue capture



Magnum Style
UMRP-1820 18g x 20cm
UMRP-1825 18g x 25cm

- Etched Cutting notch for better tissue retention
- Echo enhanced tip for ultrasound placement
- Spacer clip for easy loading

Endoscopic Urology Graspers



- Oval cup with needle
- Standard round cup
- Round cup with needle
- Rat tooth retriever
- Long oval cup
- serrated cup with needle
- Alligator jay retriever
- Serrated cup

Urovance, LLC. repairs and sells quality remanufactured endoscopes and is not associated or affiliated in any way with Olympus®, Karl Storz® or Pentax® "All product names and images represented here are trademarks or registered trademarks of their respective holders."

Contact Ben Cooley to discuss your repair or purchase needs
214.505.3896 info@urovancemedical.com