


# Silica / Silica Fibers



## Free-Beam Fiber Delivery System Specification Worksheet

<i>Holmium-ND:YAG-Diode</i>	<b>272</b>	<b>365</b>	<b>550</b>	<b>940</b>
Single-Use SMA905 Connector	UMSU-272	UMSU-365	UMSU-550	UMSU-940
Reusable SMA905 Connector	UMRU-272	UMRU-365	UMRU-550	UMRU-940
Single-Use Trimedyne™ Connector	UMTS-272	UMTS-365	UMTS-550	UMTS-940
Reusable Trimedyne™ Connector	UMTR-272	UMTR-365	UMTR-550	R UMTR B-940
Outer Diameter	400µm	730µm	750µm	1400µm
Core	272µm	365µm	550µm	940µm
Working Channel	0.5mm (1.5 Fr.) Required	1.0mm (3.0 Fr.) Required	1.0mm (3.0 Fr.) Required	1.6mm (5.0 Fr.) Required

## Surgical Fibers Quality Accessories

<i>Model Number</i>	<i>Description</i>
FS-272	Stripper for 272µm Core Fiber
FS-365	Stripper for 365µm Core Fiber
FS-550	Stripper for 550µm Core Fiber
FS-940	Stripper for 940µm Core Fiber
FC-001	Fiber Cleaver

*Single use and reusable laser fibers*